




Problem Set 11 - Partial Answers

Problem R-310 ($C_{11}H_{10}Cl_2F_2$). Interpret the 56.4 MHz ^{19}F NMR spectrum below (CCl_4 solvent). Determine the chemical shifts of the fluorines, and estimate the various coupling constants. Consider conformations of the cyclobutane ring (*J. Am. Chem. Soc.* **1962**, *84*, 2935).




Problem R-311 ($\text{C}_6\text{H}_3\text{D}_3$). Assign the peaks in the ^{13}C NMR spectrum below. The spectrum is not ^1H decoupled. Estimate the coupling constants (*J. Am. Chem. Soc.* **1967**, 88, 2967).



ANSWER

There are four kinds of carbons in this molecule




1: The only significant coupling is $^1J_{\text{CD}} = 23 \text{ Hz}$. The $^2J_{\text{CH}}$ will be 1 Hz, too small to resolve




2: This carbon will be a double intensity dd, with $^1J_{\text{CH}} = 159 \text{ Hz}$ and $^3J_{\text{CH}} = 7 \text{ Hz}$. There will also be a $^2J_{\text{CD}}$ of about 1 Hz, too small to resolve



3: Carbon 3 will be a double intensity 1:1:1 triplet of 1:1:1 triplets, $^1J_{\text{CD}} = 23 \text{ Hz}$ and $^3J_{\text{CD}} = 1 \text{ Hz}$





4: Carbon 4 will be a 1:1:1 triplet of 1:2:1 triplets, $^1J_{\text{CD}} = 23 \text{ Hz}$ and $^3J_{\text{CH}} = 7 \text{ Hz}$ the outer lines of the triplets are too small to be visible



These outer triplet peaks will be approximately 1/16 of the intensity of the superimposed central peaks, hence not detectable at this signal to noise

Problem R-28C ($C_{21}H_{32}O_3SSi$). Use the 500 MHz ($CDCl_3$) homonuclear decoupled spectra below to assign the protons of the compound shown. Source: Mark Matulenko/Burke



Problem R-28D (C₁₃H₁₄O₂)
 500 MHz ¹H NMR in CDCl₃
 Source: Wilds/C. Fry (C82)

Assign all protons in this spectrum, using the 500 MHz ¹H NMR spectrum, and the 300 MHz COSY spectrum. The 300 MHz ¹H spectrum is also provided. Explain specifically why some of the peaks are more complicated in the 300 MHz compared to the 500 MHz spectrum. Draw a conformation, and label with chemical shifts.

