


Problem R-98B ($C_{10}H_{12}O_2$). Determine the structure (or part structure) of R-98B from the ¹H NMR, ¹³C NMR and IR spectra provided.

- (a) DBE___
- (b) What information can you obtain from the IR spectrum?


(c) Interpret the 13 C NMR spectrum. The DEPT 135 spectrum shows all CH and CH $_3$ peaks as positive, and CH $_2$ peaks negative, the DEPT 90 is CH only. Identify what kind of carbon each signal corresponds to, and write possible part structures.

-	-	
No	ppm	Type of C (e.g. sp ³ CH ₂) and/or part structures (e.g. N-CH ₂)
1	146.3	
2	143.6	
3	137.6	
4	131.5	
5	120.8	
6	115.0	
7	114.2	
8	111.0	
9	55.3	
10	39.4	

(d) Analyze the multiplet at δ 5.9 in the 270 MHz ¹H NMR spectrum (shown below). Report multiplicity, coupling constants and part structure you could obtain from the signal. Draw a coupling tree.




(e) The peaks between δ 6 and δ 7 have been provided both in acetone-d6 and in benzene-d6 solution. Do an analysis of the pattern in each solvent (give J and δ), and draw the part structure identified by the pattern. Label the part structure to indicate which signals correspond to each proton. (Do not attempt to explain the solvent effect.)





Part structure:


(f) Analyze the remainder of the signals in the 1H spectrum. Label structure fragments with J and δ values.




(g) The data provided do not allow a complete structure assignment. Draw all likely structures below. Do a chemical shift calculation for the proton you have assigned to the multiplet at δ 6.6 in the acetone spectrum for each possible structure, and label the structure with the calculated value (do not try to calculate all of the proton shifts). Circle the most likely structure.









Problem R-98B (C₁₀H₁₂O₂). Determine the structure (or part structure) of R-98B from the ¹H NMR, ¹³C NMR and IR spectra provided.

(a) DBE 5

2

3

(b) What information can you obtain from the IR spectrum?

3600 cm⁻¹ sharp OH - perhaps intramolecularly H-bonded

No carbonyl 1600-1800 cm⁻¹

No acetylene, allene

(c) Interpret the ¹³C NMR spectrum. The DEPT 135 spectrum shows all CH and CH₃ peaks as positive, and CH₂ peaks negative, the DEPT 90 is CH only. Identify what kind of carbon each signal corresponds to, and write possible part structures.

Type of C (e.g. sp³ CH₂) and/or part structures (e.g. N-CH₂) No ppm

146.3 1 2 143.6 Arom, vinyl CH 3 137.6 131.5 4 Arom, vinyl CH 5 120.8 6 115.0 Arom, vinyl CH 114.2 7 Arom, vinyl CH 111.0 8

There are 11 C-H identified from the multiplicities of the ¹³C NMR spectrum. There is one H missing - must be an O-H (also seen in IR)

There are 10 different carbons - so no symmetry

- (d) Analyze the multiplet at δ 5.9 in the 270 MHz ¹H NMR spectrum (shown below). Report multiplicity, coupling constants and part structure you could obtain from the signal. Draw a coupling tree.
 - One vinyl proton (from δ)


CH₃-O


C-CH₂-C

- ddt, J = 17, 10, 7 Hz

this defines:








9

10


55.3

39.4

(e) The peaks between δ 6 and δ 7 have been provided both in acetone-d6 and in benzene-d6 solution. Do an analysis of the pattern in each solvent (give J and δ), and draw the part structure identified by the pattern. Label the part structure to indicate which signals correspond to each proton. (Do not attempt to explain the solvent effect.)




 \cdot The sunstituents must be electrom donating, from δ




MeO

(f) Analyze the remainder of the signals in the ^{1}H spectrum. Label structure fragments with J and δ values.



(g) The data provided do not allow a complete structure assignment. Draw all likely structures below. Do a chemical shift calculation for the proton you have assigned to the multiplet at δ 6.6 in the acetone spectrum for each possible structure, and label the structure with the calculated value (do not try to calculate all of the proton shifts). Circle the most likely structure.



HC

7

5

